

The top view shows a rectangular bathtub with overall dimensions of 60" in length and 36" in width. The interior length is 50" and the interior width is 28". The depth of the tub is 19 1/2". The tub features a central drain (FJ) and two overflow drains (FJ) on the left side. The tub is equipped with a water pump (WP PUMP) and a main pump (MS PUMP) on the right side. The tub is shown with a 10" wide overflow area on the left and a 2" wide overflow area on the right. The tub is shown with a 6 1/2" wide overflow area on the left and a 2" wide overflow area on the right.



The side view shows the tub's profile with a 4" wide overflow area on the left and a 2" wide overflow area on the right. The tub has a 15 1/2" high overflow area on the left and a 18 1/2" high overflow area on the right. The tub is shown with a 43 3/16" wide base and a 20" high backrest. The tub is shown with a 120° angle on the backrest and a 95° angle on the side panel. The tub is shown with a 15 1/2" high overflow area on the left and a 18 1/2" high overflow area on the right.

Note: Equipment may extend beyond the outside perimeter of the bath. Consult factory for details.

Note: Dimensional Tolerances: +1/4"

Drop-in Hole Size:	58" x 34"
Operating Capacity:	37 Gal.
Overflow Capacity:	56 Gal.
Floor Loading / sq.ft.:	49.2 lbs.
Carton Wt.:	132 lbs.
Cube:	31.3 ft. ³

We provide products listed with Intertek ETL SEMKO (Edison Testing Laboratories) U.S. and Canadian and with International Association of Plumbing and Mechanical Officials (IAPMO). Consult factory with your listing requirements.

Standard	 <p>A 120V system that includes a MicroSilk® unit; a quiet 1-1/2 HP, 1-speed whirlpool pump; an electronic control; and a 20-minute timer. The keypad also operates the optional Jetted Neck Pillow.</p>
Premier	 <p>A 120V system that includes a MicroSilk® unit; a quiet 1-1/2 HP, 1-speed whirlpool pump; bath heater; (2) Chromatherapy mood mighths; Automatic Ozone System; a digital control with LCD display; enhanced user control options; and a 20-minute timer. Control keypad operates the MicroSilk® and Whirlpool pumps, LED Chromatherapy mood lights, optional Jetted Neck Pillow, and shows equipment status, and time remaining.</p>